RGF

aquamax

TIEDRA

iLents Care, S.L.
Pintores, 5
28923 Alcorcón
MADRID ESPAÑA
www.tiedra.net

TIEDRA CE 0318

STERILE A

27380-4

ESPAÑOL - ES

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Texto revisado: 23/05/2024

INGLÉS - EN

- Read the instructions carefully and keep the leaflet for later reference

Aquamax® RGF, by TIEDRA®

How to use: To clean and preserve your RGF contact lenses, follow the next steps:

Step 1. Before handling your lenses, Wash and rinse your hands well.

Step 2. Shake before use. Apply 2-3 drops of Aquamax® RGF to both sides of the lens. Gently rub it with your thumb and forefinger for 15 seconds.

Step 3. Rinse the lens well with preservative solution or saline solution. You can place the lens in the eye or store it in the case immersed in Aquamax RGF Moist.

Step 4. In the morning, you can put your lenses on directly.

Composition: Sodium Chloride, Lauryl Ether Sulfate, Silica Gel, Titanium Dioxide, Boric Acid, Sodium Tetraborate, Disodium Edetate 0.1%, y Polyhexamethylene Biguanide 0.001% in Purified Water.

How does Aquamax® RGF work?: Aquamax® RGF is formulated with two main ingredients. The first of these is a cleaning agent capable of loosening deposits that stick to your contact lens when you rub the lenses. The second ingredient captures the loose deposits. Finally when you rinse the lenses with sterile saline or Aquamax RGF Moist, the cleaning agent and all captured deposits are removed.

What Benefits does Aquamax® RGF provide?: Aquamax® RGF performs a deep cleaning of the surface of your rigid permeable to gases contact lens, eliminating deposits of substances that adhere during use. In this way, it prevents the growth of microorganisms due to the presence of deposits, as well as the appearance of allergies to adhering substances, so the use of your contact lenses will not be shortened by this type of circumstance.

Indications (Uses): Aquamax® RGF is an aqueous solution for RGF contact lenses cleaning.

Contraindications (When do not use): Do not use with soft contact lenses. Do not use if you are allergic to any of its

components. Do not use the product if the packaging is not intact, or the cap guarantee seal is open. In this situation, contact the health center where you purchased this product.

WARNINGS: 1. Do not put Aquamax® RGF in direct contact with the eyes. After cleaning your contact lenses with Aquamax® RGF, you must rinse the lenses with Aquamax® RGF Moist product or sterile saline solution before putting them in the case immersed in Aquamax® RGF Moist.

2. PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS MAY CAUSE CORNEAL INFECTIONS OR ULCERS AND CAUSE VISION LOSS. Failure to dispose of the solution in the case after each use or leaving the lenses in tap, mineral or distilled water or any non-sterile solution may contaminate the lenses and cause eye damage and possible vision loss.

Please review the following important safety information:

Important safety information: Always follow the instructions for use. See your eye care professional regularly (at least every 6 months) and follow their directions. Wash and dry your hands thoroughly before handling lenses. Do not use tap water, mineral water or saliva to clean the lenses or to fill the case. Use only Aquamax® RGF Cleaner from the bottle to clean the lenses. Always throw away any leftover solution from the kit after each use. Do not add Solution to lens case with used solution. Clean, rinse and air dry the lens case after each use. Turn the case upside down to drain off excess solution as it dries. Replace the lens case with a new one frequently, at least when starting a new product pack. Other products such as saline solution or rewetting drops do not disinfect lenses. Always replace the maintenance solution and lenses with new ones when appropriate.

Precautions: Avoid contacting the tip of the bottle with other surfaces in order not to contaminate the lenses. Redose canister after each use. Keep out of reach of children. Consume within 90 days, once the container is opened. Keep the product with the bottle tightly closed, in a dry place, at temperatures not exceeding 30°C. You must not use the product after the expiration date that appears on the box and bottle.

Adverse reactions (problems and what to do): You may also experience the following reactions: stinging, burning or itching (irritation), reduced comfort compared to the first time you put the lenses on, a feeling of something in the eye (a foreign body, a scratch in a certain area), tearing eye wear, unusual discharge from the eye, red eyes, decreased visual acuity (reduced vision), blurred vision, glare or halos around objects, sensitivity to light (photophobia), or dry eyes.

If you experience any of these reactions: **Remove your**

lenses immediately and if the discomfort or problem goes away, carefully examine the lenses up close. If the lens is damaged, do not put it back on your eye. Put the lens back in its case and contact your eye care professional. If the lens is dirty, has chips or a foreign body, or if the problem goes away and the lens is undamaged and has been thoroughly cleaned, rinsed and disinfected, put the lens back on. **If the problem persists, remove your lenses and see your eye care professional immediately.**

Adverse effects notification: You can notify us of medical incidents or adverse reactions at the web address: web: <http://ae.tiedra.net/> You should also notify to the Health Authorities at your country.

Throw away: Dispose of the product in the usual way. Try to recycle cardboard carton, plastic and aluminium packaging as indicated in your city

Formats: Aquamax® RGP is presented in different formats. The batch number and expiration date are listed on both the product box and canister and on the neutralizing tablets.

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