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Texto revisado: 23/05/2024

OPERATING AND MAINTENANCE INSTRUCTIONS

IMPORTANT INFORMATION FOR USERS OF RIGID GAS-PERMEABLE CONTACT LENSES

CE 0318	 0318
Medical device (5.7.7)	
Do not use if packaging is damaged	
GREEN POINT	
MANUFACTURER	
LOT	
DATE OF MANUFACTURE	
EXPIRY DATE	

 iLents Care, S.L
Pintores, 5
28923 Alcorcón
MADRID ESPAÑA
www.tiedra.net

READ THE INSTRUCTIONS CAREFULLY AND KEEP THE PACKAGE LEAFLET FOR FUTURE REFERENCE.

TIEDRA® rigid gas-permeable contact lenses

- These contact lenses have been recommended by your eye care professional, it is essential that you follow their recommendations at all times.
- He will indicate the most suitable replacement time and number of hours of use.
- Visit your eye care professional regularly (at least every 6 months) and follow their instructions.
- Never wear contact lenses that have not been recommended by a professional.
- If in doubt, consult your eye care professional.
- If these contact lenses are to be used by a minor, the mother and/or father must pay special attention to all the information in these instructions for use to ensure the safe use of this contact lens by their child.
- Always follow these instructions for use.

INSTRUCTIONS FOR USE FOR RIGID GAS PERMEABLE CONTACT LENSES

EN

TIEDRA

HOW TO WEAR:

Follow the steps below to put in and take out your rigid gas permeable contact lenses.

TO PUT IN CONTACT LENSES

1. Alexa rigid gas permeable contact lenses have been supplied to your lens adapter in tamper-evident packaging. Check the condition of new packaging for signs of damage. Do not use the lenses if the original packaging is broken, damaged or opened.
2. Wash and rinse your hands thoroughly and dry them with a lint-free towel.
3. Remove the lens from the case with the tip of your index finger. Always proceed in the same order.
4. Immediately before insertion into the eye, rinse the inner side of the lens with maintenance solution and moisten it with a little moisturiser.
5. Place the lens on a lightly moistened index fingertip.
6. Be careful. The best system is:
right lens - left hand. // left lens - right hand.
7. First fit the lens that corresponds to the eye with the greater visual defect.
8. During the adaptation period, you place the lens in front of a mirror placed horizontally on the table. Good lighting is important.
9. Bring the left hand over the forehead, holding the upper eyelid with one finger. The eye is wide open. The iris is exposed.
10. Carefully place the lens in the centre of the eye, pulling the lower eyelid slightly downwards with the middle finger of the right hand.
11. Check this operation in the mirror. Carefully release your eyelids and close your eye briefly. If you now have good vision, you have succeeded. You can forget about your contact lenses for the rest of the day.

At first, you will not be able to avoid a flickering effect when inserting the lens. For this reason, you should hold your eyelids firmly in place.

If you have not positioned the lens exactly in the centre of the eye, you centre it carefully with your fingers and eyelids. If the lens falls off the finger when inserting it, it should be cleaned again by rinsing and moistening it again.

TO REMOVE CONTACT LENSES

Not all eyes are the same. For this reason, different methods are described here.

What is important is that the lens is normally centred, being approximately in the centre of the eye.

1st METHOD.

1. Tilt the head over a mirror.
2. The right lens is removed with the right hand, the left lens with the left hand.
3. Open the eyes wide so that the lens is no longer held by the eyelids.
4. Place the tip of the finger on the outer corner of the eyelid. Be careful! The fingernail must not touch the eye.
5. Hold the eye open and pull both eyelid rims taut by pulling them towards the ear. The lens falls out of the eye. Clean it immediately and put it back in the case - be careful, danger of confusing them!

2nd METHOD.

1. Check in the mirror.
2. One finger holds the upper eyelid, one finger of the other hand holds the lower eyelid in the vicinity of the outer corner of the eyelid.
3. With the help of the fingers, the eye is opened wider. The lens is free of the eyelids.
4. Then, the eyelids are simultaneously tightened in the direction of the ear. While you direct your gaze slightly towards your nose.
5. Due to the pressure of the eyelids on the edge of the lens, the lens jumps out of the eye.

WARNINGS : Problems with contact lenses and lens care products may cause corneal infections or ulcers and result in vision loss. If you do not discard the lens holder solution after each use or if you leave the lenses in tap, mineral or distilled water or any non-sterile solution, the lenses may become contaminated and cause eye damage and possible vision loss. Please review the following Important Safety Instructions for the safe use of your contact lenses.

Important safety instructions. For safe use of your contact lenses, please follow the instructions below for the maintenance of these contact lenses:

- These contact lenses should be kept in clean lens holders fully immersed in an RPG lens preservative solution to be renewed daily.
- Always use the cleaning and care solutions recommended by your eye care professional for your contact lenses. These contact lenses can be cleaned with a cleaning solution and after cleaning put in lens holders in a preservative solution or alternatively with peroxide disinfection systems for contact lenses. Your eye care professional will be able to recommend other contact lens cleaning products to suit your specific needs.
- Wash and dry your hands thoroughly before handling the lenses.
- **Do not use tap water, mineral water or saliva to clean the lenses or to refill the lens holder.**
- Always discard the remaining solution from the lens holder after each use.
- Do not add maintenance solution to the lens holder with used solution.
- Clean, rinse and air dry the lens holder after each use.
- Place the lens holder upside down to drain excess solution as it dries.
- Replace the lens holder with a new one frequently, at least when starting a container of maintenance solution.
- Replace these contact lenses and with new lenses when appropriate according to the replacement period prescribed by your eye care professional.

COMPOSITION :

Trade name of the Lens you use	Composition
ALEXA 20 ALEXA AS 20	Pasifocón A; Boston® ES; Boston® IV
ALEXA 30 ALEXA AS 30 ALEXA K 30	FluoroPerm® 30; Optimum® Classic
ALEXA 60 ALEXA AS 60 ALEXA K 60 ALEXA KC ALEXA KC XL ALEXA PRESBI ALEXA MF	HDS®; Optimum® Comfort; Boston® EO
ALEXA 90 ALEXA AS 90 ALEXA ES, ALEXA ES ADVANCE, ALEXA ES PRO, ALEXA ES 12.7 ALEXA ORTO, ALEXA AR, ALEXA AR PLUS, ALEXA AR Toric, ALEXA AR PLUS Toric	HDS®100; Optimum® Extra; Boston® XO

HOW DO THESE CONTACT LENSES WORK? :

Trade name of the Lens you use	Explanation of how it works
ALEXA 20, ALEXA 30, ALEXA 60, ALEXA 90, ALEXA AS 20, ALEXA AS 30, ALEXA AS 60, ALEXA AS 90, ALEXA K 30, ALEXA K 60, ALEXA PRESBI, ALEXA MF, ALEXA ES, ALEXA ES ADVANCE, ALEXA ES PRO, ALEXA ES 12.7, ALEXA KC, ALEXA KC XL	These contact lenses are manufactured with a geometry that is defined by the anatomical measurements of your ocular surface and the physiological requirements of your eye. With this information specific to your eye, we generate a geometric shape in the contact lens capable of modifying the direction of the light rays, so that they focus on your retina, and thus your vision is clear.
ALEXA ORTO ALEXA AR, ALEXA AR PLUS, ALEXA AR Toric, ALEXA AR PLUS Toric, ALEXA AR PRO, ALEXA AR PRO Plus, ALEXA AR H	The mechanism of action of the ALEXA ORTO lens is nocturnal orthokeratology, whose objective is to vary the curvature of the cornea to reduce its refractive power, for which the adapted lens must have certain characteristics in its design, these lenses are called inverse geometry lenses.

WHAT ARE THE BENEFITS OF CONTACT LENSES? :

Trade name of the Lens you use	Benefits you get by using these contact lenses
ALEXA 20, ALEXA 30, ALEXA 60, ALEXA 90, ALEXA AS 20, ALEXA AS 30, ALEXA AS 60, ALEXA AS 90, ALEXA K 30, ALEXA K 60, ALEXA PRESBI, ALEXA MF, ALEXA ES, ALEXA ES ADVANCE, ALEXA ES PRO, ALEXA ES 12.7, ALEXA KC, ALEXA KC XL	By using these contact lenses, your ametropic vision will be compensated, so that you will achieve a visual acuity (clear vision) that allows you to carry out your daily activities normally.
ALEXA ORTO ALEXA AR, ALEXA AR PLUS, ALEXA AR Toric, ALEXA AR PLUS Toric, ALEXA AR PRO, ALEXA AR PRO Plus, ALEXA AR H	By wearing these lenses at night, your ametropia of vision will be compensated, so that you will achieve visual acuity (clear vision) during the day, allowing you to carry out your daily activities normally.

INTENDED PURPOSE :

Trade name of the Lens you use	Intended Purpose
ALEXA 20, ALEXA 30, ALEXA 60, ALEXA 90 ALEXA AS 20, ALEXA AS 30, ALEXA AS 60, ALEXA AS 90, ALEXA K 30, ALEXA K 60, ALEXA PRESBI, ALEXA MF, ALEXA ES, ALEXA ES ADVANCE, ALEXA ES PRO, ALEXA ES 12.7, ALEXA KC, ALEXA KC XL, ALEXA ORTO, ALEXA AR, ALEXA AR PLUS, ALEXA AR Toric, ALEXA AR PLUS Toric, ALEXA AR PRO, ALEXA AR PRO Plus, ALEXA AR H	Compensation of ametropia of vision.

INDICATIONS (USES) :

Trade name of the Lens you use	Indication
ALEXA 20 ALEXA 30 ALEXA 60 ALEXA 90 ALEXA AS 20 ALEXA AS 30 ALEXA AS 60 ALEXA AS 90	Myopia, hyperopia and astigmatism compensation in adults and children with healthy corneas.
ALEXA PRESBI, ALEXA MF	Presbyopia compensation in adults.
ALEXA K 30, ALEXA K 60 ALEXA ES, ALEXA ES ADVANCE, ALEXA ES PRO, ALEXA ES 12.7, ALEXA KC, ALEXA KC XL	Adults with cornea with keratoconus. Post keratoplasty cornea, post refractive surgery cornea.
ALEXA ORTO, ALEXA AR, ALEXA AR PLUS, ALEXA AR Toric, ALEXA AR PLUS Toric, ALEXA AR PRO, ALEXA AR PRO Plus, ALEXA AR H	Myopia compensation in children and adolescents.

CONTRAINDICATIONS (REASONS NOT TO USE THIS PRODUCT): You should not use contact lenses in the following situations. If in doubt, consult your optician or ophthalmologist:

1. When ocular pathologies exist, such as (but not limited to): Ocular infections, Conjunctivitis or inflammation of the conjunctiva, Keratitis or corneal inflammation, Dry Eye, Acute inflammation of the anterior segment of the eye (iritis), Reduced corneal sensitivity, Ptosis, Chronic ocular allergy, red eye or ocular irritation.
2. Allergic reactions on the ocular surface, eyelids or eyelashes that may be caused or exaggerated by the use of contact lenses.
3. Systemic pathologies that may affect the eye or that may be exaggerated by the use of contact lenses.
4. Some medications, especially eye medications, are contraindicated in contact lens wear. If this is your situation, consult your optician or ophthalmologist.
5. Wearers who, for whatever reason, are unable to maintain or handle contact lenses properly.

PRECAUTIONS AND THINGS TO CONSIDER WHEN WEARING CONTACT LENSES:

- These contact lenses have been recommended by your eye care professional, it is essential that you follow his recommendations at all times. He or she will advise you on the most appropriate replacement time and number of hours of wear.
- Visit your eye care professional regularly (at least every 6 months) and follow their instructions.
- Never wear contact lenses that have not been recommended by a professional.
- If in doubt, consult your eye care professional.
- If these contact lenses are to be used by a minor, the mother and/or father must pay special attention to all the information in these instructions for use to ensure the safe use of this contact lens by their child.
- Always follow these instructions for use.
- Do not wear contact lenses longer than your eye care professional recommends.
- Do not use the lenses if the original vial is broken, damaged or opened. In this case, contact the contact lens fitting centre.
- The original outer packaging should be kept until the last lens of its contents has been used.
- If you do not wear your contact lenses for a few days, ask your optician for advice on proper care.
- Do not practise water sports with lenses on without wearing waterproof goggles.
- Put in contact lenses before applying eye make-up and remove them before removing make-up. Use cosmetic products recommended for contact lens wearers. Make-up can discolour your lenses.
- Avoid contact with soaps, creams, shampoos, hairsprays, or other cosmetics as they may irritate your eyes. Close your eyes if you use hairspray, as it can damage your lenses.

- Be careful when handling the lenses, excessive pressure can deform the lenses.
- Scratches affect the wearing comfort. For this reason, do not slide the lenses on hard surfaces. The best way to lift the contact lens is to touch it lightly with a moistened fingertip.
- **If you wear contact lenses during the day, DO NOT sleep with them.** Remember to remove your contact lenses before going to sleep.
- **If you wear orthokeratology lenses** to sleep with them, **NEVER use them to see with your eyes open.** Orthokeratology lenses are designed to be worn with your eyes closed while you sleep. If you wear them with your eyes open, you can damage your eyes.
- Keep contact lenses and care products out of the reach of children.
- Do not use the product after the expiry date stated on the carton and bottle.
- If your contact lenses come into contact with harmful vapours, chemicals or dangerous environments while you are wearing them, you should remove them immediately.

ADVERSE REACTIONS (PROBLEMS AND WHAT TO DO) : While wearing contact lenses and lens care solutions, you may experience the following reactions: stinging, burning or itching (irritation), reduced comfort compared to when you first wore the lenses, sensation of having something in your eye (a foreign body, a scratch in a certain area), excessive tearing of the eye, unusual discharge from the eye, red eyes, corneal oedema, corneal neovascularisation, abnormalities of the tarsal conjunctiva, iritis, conjunctivitis, reduced visual acuity (reduced vision), blurred vision, flashes or halos around objects, sensitivity to light (photophobia) or dry eyes.

If you experience any of these reactions: **Remove the lenses immediately and**

- If the discomfort or problem disappears, examine the lenses closely.
- If the lens is damaged, do not put it back in your eye. Return the lens to its case and contact your eye care professional.
- If the lens is dirty, has an eyelash or foreign body, or if the problem disappears and the lens is undamaged and you have thoroughly cleaned, rinsed and disinfected it, put the lens back in.
- **If the problem persists, remove the lenses and see your eye care professional immediately.**

ADVERSE EVENT REPORTING:

You can report medical incidents or adverse reactions to us at the following web address:

<http://ae.tiedra.net/>.

You can also report to the Health Authorities in the country where you are located.

Disposal: Dispose of the product in the normal way. Try to recycle cardboard and plastic as indicated in your city.

Presentation: the product is available in different formats. The batch number and expiry date can be found on the product container.

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