ESPAÑOL - ES

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**Eye drops
in single-dose containers**
Hypertonic, sterile, ophthalmic solution

Contiene **25** envases
unidosis de 0,4 ml

Contains **25** single-dose
containers of 0.4 ml

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Texto revisado: 25/09/2025

INGLÉS - EN

- Read the instructions carefully and keep the leaflet for later reference

OCUMAX 5®, by TIEDRA®

How does hypertonic sodium chloride (or brand name) work?

Hypertonic sodium chloride (or brand name) is a hypertonic sodium chloride solution that acts as a semi-permeable hypertonic membrane when it comes into contact with the corneal epithelium. Corneal oedema is characterised by excess fluid in the cornea. Hypertonic sodium chloride (or brand name) extracts excess fluid from the cornea through osmosis, helping to keep the cornea clear and free from inflammation, thereby improving vision quality.

Indications

Hypertonic sodium chloride (or brand name) is indicated for the reduction of corneal oedema, which can be caused by various factors such as bullous keratitis, post-cataract surgery and hereditary corneal dystrophy or Fuchs' dystrophy.

Benefits

CSH will reduce corneal oedema by removing excess fluid, thereby improving transparency and consequently improving your visual quality.

Contraindications

Do not use Hypertonic Sodium Chloride (brand name/trade name) in the following cases:

- **Hypersensitivity or known allergy to any of the components of the solution.**
- **Open wounds or severe epithelial lesions of the cornea**, unless expressly indicated by a doctor.
- If in doubt, always consult your doctor or pharmacist before using the product.

How to use

Hypertonic sodium chloride (or brand name) is used ophthalmically (applied to the eye).

The recommended dose is one or two drops in the affected eye or eyes every six to eight hours. The number of daily applications and the duration of treatment may be modified according to medical criteria.

The solution should be used immediately after opening.

Wash your hands before use.

For correct administration, follow these instructions:

1. Remove a single-dose container from the strip.
2. Check that the single-dose container is not damaged.
3. Hold the single-dose container between the index finger and thumb of one hand.
4. Using the index finger and thumb of your other hand, gently lever the winged piece to open it. Do not twist to break.
5. Remove the winged piece.
6. Gently pull down the lower eyelid to form a pocket. Turn the single-dose container upside down and squeeze it to apply one or two drops into the eye. Keep your eye closed for one (1) or two (2) minutes.

Do not reuse the single-dose container even if there is solution remaining, as the remaining solution could become contaminated and cause infection.

If you are using more than one ophthalmic product, consult your doctor or pharmacist about how long to wait between applications of each product.

WARNINGS

Warning: this product is designed for use in eyes with corneal oedema. Do not use this product in an eye without oedema, as it could cause potentially dangerous damage to the surface of the eye.

PRECAUTIONS

- Do not use if the single-dose container is broken or damaged, as the product may be spoiled or contaminated.
- For ophthalmic use only, do not ingest.
- Avoid contact between the tip of the dropper and the surface of the eye or eyelid to prevent possible eye irritation from contact between the tip of the dropper and these surfaces.
- Avoid contact between the tip of the dropper and the surface of the eye or eyelid, as well as fingers, hands or non-sterile tissues to prevent contamination of the ophthalmic solution.
- Keep the product out of sunlight and at a room temperature below 25°C.
- Discard the opened single-dose container after use, as the remaining solution could become contaminated and its use could cause infection.
- Do not use after the expiry date indicated on the box and on each single dose.
- In some cases, a burning sensation or irritation may be experienced after instillation, which disappears within a few minutes. If this persists, stop using the product and consult your doctor or pharmacist.
- Keep out of the sight and reach of children.

Adverse reactions

This product is a sterile hypertonic solution and, as such, its use may cause adverse effects such as

Common adverse effects:

- Transient stinging or burning sensation immediately after instillation.
- Mild eye irritation (itching, redness, tearing).
- Foreign body sensation or slight eye discomfort.

Less common adverse effects:

- Eye pain (if irritation is severe or prolonged).
- Temporary blurred vision after application (due to osmotic changes and transient increase in tearing).

Rare adverse effects:

- Mild keratoepitheliitis or worsening of epithelial microlesions if the cornea is severely damaged.
- Increased sensation of dry eyes in predisposed patients.

COMPOSITION

Sodium chloride 5%, hypromellose, boric acid, borax and purified water.

Adverse effects. Based on our experience with OCUMAX 5®, no adverse effects are expected. If you experience any adverse effects, consult your eye care professional.

Adverse effects notification: You can notify us of medical incidents or adverse reactions at the web address: web: <https://ae.tiedra.net/> You should also notify to the Health Authorities at your country.

Throw away: Dispose of the product in the usual way. Try to recycle cardboard cartoning, plastic and aluminium packaging as indicated in your city.

Formats: The product is presented in different formats. The batch number and expiration date are listed on both the product box and label

Reviewed Text: 25/09/2025