ESPAÑOL - ES

AQUAMAX BIO

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AQUAMAX BIO

Moisturising drops with sodium hyaluronate

COMPOSITION: Sodium chloride, boric acid, sodium tetraborate, disodium edetate, hydroxyethylcellulose, hyaluronic acid, polyhexamethylene biguanide (PHMB) and purified water.

HOW IT WORKS: Aquamax® Bio combines sodium hyaluronate and hydroxyethylcellulose, two moisturising agents that form a protective film on the surface of the eye, helping to retain moisture and prolong the eye's natural lubrication.

BENEFITS: Aquamax® Bio provides rapid and long-lasting relief from dry eyes, improves comfort in situations of irritation and helps keep eyes fresh and hydrated. Its formulation is specially designed for contact lens wearers and for those who spend long periods in front of screens or in unfavourable environments.

INDICATIONS: Aquamax® BIO improves comfort and relieves the sensation of dry eyes. For cases where extra eye moisturisation is needed, such as:

- Stuffy environments with a lack of oxygen
- Dry atmosphere, wind, air conditioning
- Use of contact lenses
- Lack of sleep, stress
- Activities that strain the eyes, such as: use of screens, long journeys, prolonged reading, etc.

CONTRAINDICATIONS:

- Do not use if you are allergic to any of its components.
- Do not use if the container is not intact or the seal on the cap is broken. Contact the health centre where you purchased this product.

HOW TO USE:

Insert one or two drops directly into the eye on the inside of the conjunctiva, close your eyelids and move your eyes for a few seconds.

PRECAUTIONS:

- Do not touch the dropper with other surfaces, as this could contaminate it and cause infection.
- Keep the bottle tightly closed in a dry place away from light.

- Store at temperatures between 5°C and 30°C.
- Check that the product seal is intact before first use. Use within 60 days of opening the container.
- In case of persistent irritation, discontinue use and consult a professional.
- Do not use the product after the expiry date stated on the bottle.
- If you are using more than one ophthalmic product, consult your eye care professional about how long to wait between applications of each product.
- Keep out of the reach of children.

ADVERSE EFFECTS: Based on our experience with AQUAMAX® Bio, no adverse effects are expected. If you experience any adverse effects, consult your eye care professional.

ADVERSE EFFECTS NOTIFICATION: You can notify us of medical incidents or adverse reactions at the web address: web: <https://ae.tiedra.net/>

You should also notify to the Health Authorities at your country.

THROW AWAY: Dispose of the product in the usual way. Try to recycle cardboard cartoning, plastic and aluminium packaging as indicated in your city.

FORMATS: Aquamax® BIO is presented in different formats. The batch number and expiration date are listed on both the product box and label.

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